



# SIMPLE POINT OF CARE HbA1c and ACR

**SIMPLE**

**FAST**

**RELIABLE**

The Afinion™ AS100 Analyzer System is designed for point of care testing and is for *in vitro* diagnostic use only.

Refer to operator's manual for warnings, precautions and limitations regarding proper use of this instrument. Manuals can be found on [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com).

*afinion*™



*Put science on your side.*

 **Abbott**  
A Promise for Life



## ADA Guidelines for HbA1c

The American Diabetes Association (ADA) recommends that someone with diabetes be tested for HbA1c a minimum of two and up to four times per year. Diabetes Control and Complications Trial (DCCT) findings show that people with diabetes who keep their HbA1c levels close to 7% or below have a much better chance of preventing and/or delaying diabetic problems (i.e., heart, nerve, eye and kidney damage).

## ADA & NKF Guidelines for ACR

In addition the ADA also recommends annual screening of ACR in patients with Type 1 diabetes for a duration of  $\geq 5$  years and all Type 2 diabetes patients, starting at diagnosis and during pregnancy. The National Kidney Foundation (NKF) recommends that screening of diabetes should include measurement of Urinary ACR in a spot urine sample.

# 3-Minute HbA1c Test CLIA-Waived

## As Easy as 1 - 2 - 3

### 1 - DRAW



### SIMPLE

- Factory Calibrated
- Ready-To-Use Controls
- Operator-Free System Check

### FAST

- Three-Minute HbA1c Test Time

### 2 - INSERT



### RELIABLE

- National Glycohemoglobin Standardization Program (NGSP) Certified
- Accurate results traceable to the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC) and Diabetes Control and Complications Trial (DCCT)

### 3 - RUN



## Current Tests Available

- HbA1c – CLIA-waived
- Albumin/Creatinine Ratio (ACR)

For *in vitro* diagnostic use only. Refer to operator's manual for warnings, precautions and limitations regarding proper use of this instrument. Manuals can be found on [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com)

# Instrument Specifications

- Size: 12.6 x 6.7 x 6.7 in.
- Weight: 11 lbs.
- Input: 100-240 VAC, 50-60 Hz, 42 W



## CONTACT YOUR LOCAL AFINION REPRESENTATIVE:

Physician Sales and Service  
1-904-332-3000  
[www.pssd.com](http://www.pssd.com)

Abbott Diagnostics  
1-877-4ABBOTT (1-877-422-2688)  
[www.abbottdiagnostics.com](http://www.abbottdiagnostics.com)

Axis-Shield  
1-866-795-2947

## TECHNICAL SUPPORT

1-877-4-AFINION  
(1-877-423-4646)  
[techsupport@us.axis-shield.com](mailto:techsupport@us.axis-shield.com)  
24 hours a day, 7 days a week

Description	PSS List Number	Abbott List Number
Afinion™ AS100 Analyzer	181607	06L10
Afinion™ Cleaning Kit*	176281	08L90-50
Afinion™ Bar Code Reader*	176282	06L12
Afinion™ Printer*	176283	06L11
Afinion™ HbA1c Test Kit (15 tests)	176280	06L13-20
Afinion™ HbA1c Control Kit (low & high)	176279	06L13-10
Afinion™ ACR Test Kit (15 tests)	203019	06L14-20
Afinion™ ACR Control Kit (low & high)	203020	06L14-10

\*Optional Accessory

The Afinion™ AS100 Analyzer provides HbA1c testing for today's healthcare professionals. The Afinion™ Test System utilizes the latest technology and provides a simple, fast and reliable method by which to monitor patients' disease progression.

Refer to operator's manual for warnings, precautions and limitations regarding proper use of this instrument. Manuals can be found on [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com). Afinion™ HbA1c is an *in vitro* diagnostic test for quantitative determination of glycosylated hemoglobin (% hemoglobin A1c, HbA1c) in human whole blood.

*Put science on your side.* is a trademark of Abbott Laboratories in various jurisdictions. Afinion AS100 Analyzer is a registered trademark of Axis-Shield POC AS in the U.S.

[www.abbottdiagnostics.com](http://www.abbottdiagnostics.com)  
1-877-4ABBOTT

© 2009 Abbott Laboratories MS\_08\_16909/v2- 5 Nov Printed in USA

 **Abbott**  
Diagnostics